



## CERTIFICATE OF COMPLIANCE

### LLOYD'S REGISTER NEDERLAND B.V.

hereby declares that the certification assessment has demonstrated that

### COMPACT DRY BC

Manufactured and supplied by:  
Nissui Pharmaceutical Co. Ltd.  
20th Floor, 3-24-6 Ueno Frontier Tower  
Ueno, Taito-ku, Tokyo, 110-8736, JAPAN

has been validated and revealed to be at least equivalent to the reference method as demonstrated by the validation study report. The summary of the validation report is available on the MicroVal website:  
[www.microval.org](http://www.microval.org).

Reference method: ISO 7932:2004 Microbiology of food and animal feeding stuffs - Horizontal method for the enumeration of presumptive *Bacillus cereus* – Colony count technique at 30°C

Scope: a broad range of foods

The validation and certification has been performed in accordance with EN ISO 16140-2: 2016 and the MicroVal Rules and Certification Scheme version 8.

Limitations of the method: During the inclusivity study two strains (*Bacillus cytotoxicus* and *Bacillus pseudomycolides*) from the *Bacillus cereus* group were not detected by the alternative method whereas they were detected by the reference method.

Certificate no.: 2019LR87

First approval date: 6 March 2020

Expiry date: 5 March 2024



ISSUED BY: Lloyd's Register Nederland B.V.  
Rotterdam, The Netherlands