

○ライサスS4で同定可能なグラム陰性菌

ID	菌名コード	菌名	菌名2
1	1	Ced. davisae	Cedecea davisae
2	2	Ced. lapagei	Cedecea lapagei
3	3	Cit. amalonaticus	Citrobacter amalonaticus
4	4	Cit. freundii	Citrobacter freundii
5	5	Cit. koseri	Citrobacter koseri
6	6	Edw. tarda	Edwardsiella tarda
7	7	Ent. aerogenes	Enterobacter aerogenes
8	8	Ent. amnigenus	Enterobacter amnigenus
9	9	Ent. asburiae	Enterobacter asburiae
10	10	Ent. cancerogenus	Enterobacter cancerogenus
11	11	Ent. cloacae	Enterobacter cloacae
12	12	Ent. gergoviae	Enterobacter gergoviae
13	13	Ent. sakazakii	Enterobacter sakazakii
14	14	Esc. coli	Escherichia coli
15	15	Esc. coli (4GLC-/SOR-)	
16	16	Esc. coli (4GLC-/SOR+)	
17	17	Esc. coli (LS-)	
18	18	Esc. coli (RHA-)	
19	19	Esc. fergusonii	Escherichia fergusonii
20	20	Esc. hermannii	Escherichia hermannii
21	21	Ewi. americana	Ewingella americana
22	22	Haf. alvei	Hafnia alvei
23	23	Kleb. ornithinolytica	Klebsiella ornithinolytica
24	24	Kleb. oxytoca	Klebsiella oxytoca
25	25	Kleb. ozaenae	Klebsiella pneumoniae/ozaenae
26	26	Kleb. pneumoniae 1	Klebsiella pneumoniae/pneumoniae
27	27	Kleb. pneumoniae 2	
28	28	Kleb. rhinoscleromatis	Klebsiella pneumoniae/rhinoscleromatis
29	29	Kleb. terrigena	Klebsiella terrigena
30	30	Klu. ascorbata	Kluyvera ascorbata
31	31	Lec. adecarboxylata	Leclercia adecarboxylata
32	32	Mor. morganii	Morganella morganii
33	33	Pan. agglomerans	Pantoea agglomerans
34	34	Pro. mirabilis	Proteus mirabilis
35	35	Pro. penneri	Proteus penneri
36	36	Pro. vulgaris	Proteus vulgaris
37	37	Prov. alcalifaciens	Providencia alcalifaciens
38	38	Prov. rettgeri	Providencia rettgeri
39	39	Prov. stuartii	Providencia stuartii
40	40	Rah. aquatilis	Rahnella aquatilis
41	41	Sal. arizonae	Salmonella arizonae
42	42	Sal. gallinarum	Salmonella gallinarum
43	43	Sal. paratyphi A	Salmonella paratyphi A
44	44	Sal. typhi	Salmonella typhi
45	45	Salmonella spp.	Salmonella spp.
46	46	Ser. ficaria	Serratia ficaria
47	47	Ser. fonticola	Serratia fonticola
48	48	Ser. liquefaciens	Serratia liquefaciens
49	49	Ser. marcescens 1	Serratia marcescens
50	50	Ser. marcescens 2	
51	51	Ser. odorifera	Serratia odorifera
52	52	Ser. plymuthica	Serratia plymuthica
53	53	Ser. rubidaea	Serratia rubidaea
54	54	Shi. dysenteriae	Shigella dysenteriae
55	55	Shi. sonnei	Shigella sonnei
56	56	Shigella spp.	Shigella spp.
57	57	Yer. enterocolitica	Yersinia enterocolitica

58	58	Yer. frederiksenii	Yersinia frederiksenii
59	59	Yer. intermedia	Yersinia intermedia
60	60	Yer. kristensenii	Yersinia kristensenii
61	61	Yer. pseudotuberculosis	Yersinia pseudotuberculosis
62	62	Yok. regensburgei	Yokenella regensburgei
63	63	Aci. baumannii	Acinetobacter baumannii
64	64	Aci. haemolyticus	Acinetobacter haemolyticus
65	65	Acinetobacter spp.	Acinetobacter spp.
66	66	Alc. faecalis	Alcaligenes faecalis
67	67	Alc. xylosoxidans	Alcaligenes xylosoxidans/xylosoxidans
68	68	Bor. bronchiseptica	Bordetella bronchiseptica
69	69	Brev. diminuta	Brevundimonas diminuta
70	70	Brev. vesicularis	Brevundimonas vesicularis
71	71	Burk. cepacia	Burkholderia cepacia
72	72	Burk. mallei	Burkholderia mallei
73	73	Burk. pseudomallei	Burkholderia pseudomallei
74	74	Chr. indologenes	Chryseobacterium indologenes
75	75	Chr. meningosepticum	Chryseobacterium meningosepticum
76	76	Flav. oryzihabitans	Flavimonas oryzihabitans
77	77	Myr. odoratus	Myroides odoratus
78	78	Pse. aeruginosa	Pseudomonas aeruginosa
79	79	Pse. alcaligenes	Pseudomonas alcaligenes
80	80	Pse. fluorescens	Pseudomonas fluorescens
81	81	Pse. mendocina	Pseudomonas mendocina
82	82	Pse. putida	Pseudomonas putida
83	83	Pse. stutzeri	Pseudomonas stutzeri
84	84	Rals. pickettii	Ralstonia pickettii
85	85	She. putrefaciens	Shewanella putrefaciens
86	86	Ste. maltophilia	Stenotrophomonas maltophilia
87	87	Aer. hydrophila	Aeromonas hydrophila
88	88	Aer. sobria	Aeromonas sobria
89	89	Ple. shigelloides	Plesiomonas shigelloides
90	90	Vib. alginolyticus	Vibrio alginolyticus
91	91	Vib. cholerae	Vibrio cholerae
92	92	Vib. fluvialis	Vibrio fluvialis
93	93	Vib. mimicus	Vibrio mimicus
94	94	Vib. parahaemolyticus	Vibrio parahaemolyticus
95	95	Vib. vulnificus	Vibrio vulnificus

○ライサスS4で同定可能なグラム陽性菌

ID	菌名コード	菌名	菌名2
1	1001	Ent. avium	Enterococcus avium
2	1002	Ent. casseliflavus	Enterococcus casseliflavus
3	1003	Ent. durans	Enterococcus durans
4	1004	Ent. faecalis	Enterococcus faecalis
5	1005	Ent. faecium	Enterococcus faecium
6	1006	Ent. gallinarum	Enterococcus gallinarum
7	1007	Mic. luteus	Micrococcus luteus
8	1008	Sta. aureus	Staphylococcus aureus
9	1009	Sta. capitis	Staphylococcus capitis
10	1010	Sta. caprae	Staphylococcus caprae
11	1011	Sta. cohnii	Staphylococcus cohnii
12	1012	Sta. epidermidis	Staphylococcus epidermidis
13	1013	Sta. haemolyticus	Staphylococcus haemolyticus
14	1014	Sta. hominis	Staphylococcus hominis
15	1015	Sta. hyicus	Staphylococcus hyicus
16	1016	Sta. intermedius	Staphylococcus intermedius
17	1017	Sta. lugdunensis	Staphylococcus lugdunensis
18	1018	Sta. saprophyticus	Staphylococcus saprophyticus
19	1019	Sta. schleiferi	Staphylococcus schleiferi
20	1020	Sta. sciuri	Staphylococcus sciuri
21	1021	Sta. simulans	Staphylococcus simulans
22	1022	Sta. warneri	Staphylococcus warneri
23	1023	Sta. xylosus	Staphylococcus xylosus
24	1024	Str. agalactiae 1	Streptococcus agalactiae
25	1025	Str. anginosus	Streptococcus anginosus
26	1026	Str. bovis	Streptococcus bovis
27	1027	Str. canis	Streptococcus canis
28	1028	Str. constellatus	Streptococcus constellatus/constellatus
29	1029	Str. dysgalactiae	Streptococcus dysgalactiae/dysgalactiae
30	1030	Str. equi	Streptococcus equi/equi
31	1031	Str. equisimilis	Streptococcus dysgalactiae/equisimilis
32	1032	Str. intermedius	Streptococcus intermedius
33	1033	Str. mitis	Streptococcus mitis
34	1034	Str. mutans	Streptococcus mutans
35	1035	Str. pneumoniae	Streptococcus pneumoniae
36	1036	Str. pyogenes	Streptococcus pyogenes
37	1037	Str. salivarius	Streptococcus salivarius
38	1038	Str. sanguinis	Streptococcus sanguinis
39	1039	Str. suis	Streptococcus suis
40	1040	Str. cristatus	Streptococcus cristatus
41	1041	Str. gordonii	Streptococcus gordonii
42	1042	Str. oralis	Streptococcus oralis
43	1043	Str. parasanguinis	Streptococcus parasanguinis
44	1044	Str. pharyngis	Streptococcus constellatus/pharyngis

○ライサスS4で測定可能な薬剤感受性グループ

ID	菌種名	18時間法グループ名
1	Citro.Koseri	Enterobacteriaceae
2	Citrobacter freundii	Enterobacteriaceae
3	Citrobacter spp.	Enterobacteriaceae
4	Enterobacter aerogenes	Enterobacteriaceae
5	Enterobacter cloacae	Enterobacteriaceae
6	Enterobacter spp.	Enterobacteriaceae
7	Escherichia coli	E.coli
8	Klebsiella oxytoca	K.oxytoca
9	Klebsiella pneumoniae	K.pneumoniae
10	Klebsiella spp.	Enterobacteriaceae
11	Morganella morganii	Enterobacteriaceae
12	Proteus mirabilis	Proteus mirabilis
13	Proteus vulgaris	Enterobacteriaceae
14	Providencia spp	Enterobacteriaceae
15	Salmonella spp.	Enterobacteriaceae
16	Serratia marcescens	Enterobacteriaceae
17	Serratia spp.	Enterobacteriaceae
18	Shigella spp.	Enterobacteriaceae
19	Yersinia spp.	Enterobacteriaceae
20	Acinetobacter baumannii	Acinetobacter spp.
21	Acinetobacter spp.	Acinetobacter spp.
22	Alcaligenes spp.	ブドウ糖非発酵性グラム陰性桿菌
23	Burkholderia cepacia	Burkholderia cepacia
24	Burkholderia spp.	ブドウ糖非発酵性グラム陰性桿菌
25	Flavobacterium spp.	ブドウ糖非発酵性グラム陰性桿菌
26	Pseudomonas aeruginosa	P.aeruginosa
27	Pseudomonas spp.	P.aeruginosa
28	Sten.maltophilia	Sten.maltophilia
29	Aeromonas spp.	Vibrio group
30	Plesiomonas spp.	Vibrio group
31	Vibrio cholerae	Vibrio group
32	Vibrio spp.	Vibrio group
33	Other Enterobacteriaceae	Enterobacteriaceae
34	Other Non-Enterobacteriaceae	ブドウ糖非発酵性グラム陰性桿菌
35	Haemophilus spp.	Haemophilus spp.
36	Micrococcus spp.	CNS
37	Sta.aureus	S.aureus
38	CNS	CNS
39	Staphylococcus spp.	CNS
40	Enterococcus avium	Enterococcus
41	Enterococcus faecalis	Enterococcus
42	Enterococcus faecium	Enterococcus
43	Enterococcus casseliflavus	Enterococcus
44	Enterococcus gallinarum	Enterococcus
45	Enterococcus spp.	Enterococcus
46	Strep.pneumoniae	S.pneumoniae
47	Strep.pyogenes	beta Streptococcus
48	Streptococcus spp.	NO-ID
49	Sta.capitis	CNS
50	Sta.epidermidis	CNS
51	Sta.haemolyticus	CNS
52	Sta.saprophyticus	CNS
53	Strep.agalactiae	beta Streptococcus
54	Strep.anginosus	viridans Streptococcus
55	Strep other S.pneumoniae1	viridans Streptococcus
56	Strep other S.pneumoniae2	beta Streptococcus
57	Sta.lugdunensis	S.lugdunensis

58	NO-ID	NO-ID
59	non-fastidious(GN)	NO-ID
60	Enterobacteriaceae	Enterobacteriaceae
61	non-Enterobacteriaceae	ブドウ糖非発酵性グラム陰性桿菌
62	non-fastidious(GP)	NO-ID
63	CPS.Other than S.aureus	S.aureus
64	fastidious(GP)	NO-ID
65	beta Streptococcus spp.	beta Streptococcus
66	viridans Streptococcus spp.	viridans Streptococcus
67	bla(-)H.influenzae	Haemophilus spp.
68	bla(+ )H.influenzae	Haemophilus spp.
69	Candida spp.	Candida spp.
70	B. fragilis group	B. fragilis group
71	non B. fragilis group	non B. fragilis group
72	Moraxella catarrhalis	Moraxella catarrhalis
73	Corynebacterium spp.	Corynebacterium spp.